

Rabbit IgG Fc fragment

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

| | |
|-------------|--|
| Description | RabbitIgGfc fragment is the crystallizable fragment (Fc) of the rabbit immunoglobulin G (IgG) molecule, suitable for use as an immunolabel in various immunoassays, including ELISA. |
|-------------|--|

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481