

RI-AG03 acetate

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Keep away from moisture

Storage:

Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description

RI-AG03 acetate is the acetate form of RI-AG03. RI-AG03 is an orally active peptide inhibitor that can cross the blood-brain barrier to inhibit Tau aggregation. It effectively prevents Tau aggregation in both in vivo and in vitro models and improves associated neurodegenerative and behavioral phenotypes. RI-AG03 is useful for researching tauopathies such as Alzheimer's disease.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481