

Umedaptanib pegol sodium

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Umedaptanib pegol is a nucleic acid aptamer targeting anti-fibroblast growth factor-2 (FGF2), and it is utilized in the research of neovascular age-related macular degeneration (nAMD).
Targets(IC50)	FGFR

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481