

BCPT02

Chemical Properties

CAS No. : 3025473-38-2
Formula: C76H111N11O33
Molecular Weight: 1706.77
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	BCPT02 is a linker for synthesizing antibody-drug conjugates (ADC).
Targets(IC50)	ADC Linker

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.5859 mL	2.9295 mL	5.859 mL
5 mM	0.1172 mL	0.5859 mL	1.1718 mL
10 mM	0.0586 mL	0.293 mL	0.5859 mL
50 mM	0.0117 mL	0.0586 mL	0.1172 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481