

Procyanidin B3-3-O-gallate

Chemical Properties

CAS No. :	86631-41-6
Formula:	C37H30O16
Molecular Weight:	730.63
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

Biological Description

Description	Procyanidin B3-3-O-gallate is a naturally occurring tannin compound found in the plant <i>Sanguisorba officinalis</i> .
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.3687 mL	6.8434 mL	13.6868 mL
5 mM	0.2737 mL	1.3687 mL	2.7374 mL
10 mM	0.1369 mL	0.6843 mL	1.3687 mL
50 mM	0.0274 mL	0.1369 mL	0.2737 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481