

Chitosan oligosaccharride (MW 5000)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Chitosan oligosaccharide (MW 5000) (COS (MW 5000)) is a degradation product of chitosan, consisting of an oligosaccharide formed by linking 2-20 units of glucosamine through β -(1 \rightarrow 4) glycosidic bonds.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481