

Mouse CB11 of type II collagen

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Mouse CB11 of type II collagen (Mouse cyanogen bromide digested peptide 11) is a cyanogen bromide cleaved peptide fragment of type II collagen, useful for investigating B cell and T cell epitopes of type II collagen. This pepsin-free product is ideal for studying T cell epitope specificity via proliferation assays and can also serve as an antigen for analyzing antibody epitope specificity through ELISA.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481