

Carbomer 990

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Carbomer 990 is an acrylic polymer used as a thickening agent. It can form a high-viscosity gel, providing a stable gel matrix, and exhibits excellent resistance to heat, light, and microbial contamination. Carbomer 990 aids in drug release and topical applications, making it valuable in the pharmaceutical field.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481