

Apazunersen sodium

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Apazunersen sodium is an antisense oligonucleotide that targets and inhibits the expression of UBE3A antisense transcript (UBE3A-AS), thereby preventing the silencing of the paternal allele of the UBE3A gene and reactivating the expression of the defective protein, used in the study of Angelman syndrome.
Targets(IC50)	E1/E2/E3 Enzyme

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481