

## Ledipasvir impurity 55

## Chemical Properties

CAS No. :	1499193-60-0
Formula:	C <sub>13</sub> H <sub>6</sub> BrF <sub>2</sub> I
Molecular Weight:	407
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

## Biological Description

Description	Ledipasvir impurity 55 is a byproduct found in Ledipasvir.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.457 mL	12.285 mL	24.570 mL
5 mM	0.4914 mL	2.457 mL	4.914 mL
10 mM	0.2457 mL	1.2285 mL	2.457 mL
50 mM	0.0491 mL	0.2457 mL	0.4914 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481