

Tepyrimidine impurity 9

Chemical Properties

CAS No. :	2101241-71-6
Formula:	C ₉ H ₁₀ ClN ₃ O ₃
Molecular Weight:	243.65
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

Biological Description

Description	Tepyrimidine impurity 9 is a variant of a Tepyrimidine impurity.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.1042 mL	20.5212 mL	41.0425 mL
5 mM	0.8208 mL	4.1042 mL	8.2085 mL
10 mM	0.4104 mL	2.0521 mL	4.1042 mL
50 mM	0.0821 mL	0.4104 mL	0.8208 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481