

Dabigatran impurity 47

Chemical Properties

CAS No. :	2225986-95-6
Formula:	C34H41N7O5
Molecular Weight:	627.75
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

Biological Description

Description	Dabigatran impurity 47 is an impurity of Dabigatran.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.593 mL	7.965 mL	15.9299 mL
5 mM	0.3186 mL	1.593 mL	3.186 mL
10 mM	0.1593 mL	0.7965 mL	1.593 mL
50 mM	0.0319 mL	0.1593 mL	0.3186 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481