

## Cibigirsen

### Chemical Properties

CAS No. : 160614-28-8

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Cibigirsen is an antisense oligonucleotide designed to target the insulin-like growth factor I receptor (IGF-1R).
Targets(IC50)	IGF-1R

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481