

Dabigatran impurity 2

Chemical Properties

CAS No. :	1580491-16-2
Formula:	C32H38N8O4
Molecular Weight:	598.71
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

Biological Description

Description	Dabigatran impurity 2 is an impurity of Dabigatran.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.6703 mL	8.3513 mL	16.7026 mL
5 mM	0.3341 mL	1.6703 mL	3.3405 mL
10 mM	0.167 mL	0.8351 mL	1.6703 mL
50 mM	0.0334 mL	0.167 mL	0.3341 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481