

Tipiracil impurity 1

Chemical Properties

CAS No. :	2069937-25-1
Formula:	C ₉ H ₁₂ N ₄ O ₂
Molecular Weight:	208.22
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

Biological Description

Description	Tipiracil impurity 1 is a type of impurity associated with Tipiracil.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.8026 mL	24.0131 mL	48.0261 mL
5 mM	0.9605 mL	4.8026 mL	9.6052 mL
10 mM	0.4803 mL	2.4013 mL	4.8026 mL
50 mM	0.0961 mL	0.4803 mL	0.9605 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481