

GPI2A

Chemical Properties

CAS No. : 2857819-69-1

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description

GPI2A is an antisense oligonucleotide that complements a specific region of the HIV-1 gag gene. It significantly inhibits p55 and its cleavage products, p39/41.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481