

HBpep-K

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Keep away from moisture

Storage:

Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description

HBpep-K is a polypeptide capable of forming coacervate microdroplets. It efficiently encapsulates and delivers small peptides, proteins, and mRNA, while also protecting mRNA from RNase degradation. HBpep-K is commonly utilized in drug delivery applications.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481