

Collagen peptide Type I

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Keep away from moisture
Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	Collagenpeptide Type I is an extracellular matrix protein widely utilized in the food, cosmetics, pharmaceutical, and biomedical industries. The product consists of collagen short peptides that have been enzymatically hydrolyzed.
Targets(IC50)	MMP

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481