

Sunitinib impurity 28

Chemical Properties

CAS No. :	1186651-51-3
Formula:	C ₁₆ H ₁₄ FN ₃ O ₂
Molecular Weight:	299.31
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

Biological Description

Description	Sunitinib impurity 28 is an impurity of Sunitinib.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.341 mL	16.7051 mL	33.4102 mL
5 mM	0.6682 mL	3.341 mL	6.682 mL
10 mM	0.3341 mL	1.6705 mL	3.341 mL
50 mM	0.0668 mL	0.3341 mL	0.6682 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481