

ALN-12115

Chemical Properties

CAS No. : 3031494-36-4

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	ALN-12115 is an siRNA targeting the microtubule motor protein (KSP) and is part of ALN-VSP, a lipid nanoparticle formulation containing two siRNAs aimed at KSP and VEGF, with potential antitumor activity.
Targets(IC50)	Kinesin

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481