

Sofosbuvir impurity 3

Chemical Properties

CAS No. :	1496551-72-4
Formula:	C10H13ClN2O5
Molecular Weight:	276.67
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

Biological Description

Description	Sofosbuvirimpurity 3 is an impurity found in Sofosbuvir.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.6144 mL	18.0721 mL	36.1441 mL
5 mM	0.7229 mL	3.6144 mL	7.2288 mL
10 mM	0.3614 mL	1.8072 mL	3.6144 mL
50 mM	0.0723 mL	0.3614 mL	0.7229 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481