

Fura-2 leakage resistant AM

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage:

Keep away from direct sunlight

Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Fura-2 leakage resistant AM (Fura-PE3) is a ratiometric calcium (calcium) indicator. It exhibits calcium ion binding properties ($K_d = 146$ nM) and fluorescent characteristics (emission wavelength: 510 nm; alternative excitation wavelengths: 340 nm and 380 nm). Fura-2 leakage resistant AM is utilized for monitoring calcium ion concentration in platelets.
Targets(IC50)	Calcium Channel

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481