


CRBN ligand-734

Chemical Properties

CAS No. :	2654830-24-5
Formula:	C ₁₂ H ₁₀ FIN ₄ O ₂
Molecular Weight:	388.14
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	CRBN ligand-734 is a ligand for the E3 ubiquitin ligase cereblon (CRBN) that facilitates the recruitment of cereblon proteins. It can be linked to a target protein ligand via a linker to form a PROTAC.
Targets(IC50)	Ligands for E3 Ligase

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.5764 mL	12.8819 mL	25.7639 mL
5 mM	0.5153 mL	2.5764 mL	5.1528 mL
10 mM	0.2576 mL	1.2882 mL	2.5764 mL
50 mM	0.0515 mL	0.2576 mL	0.5153 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481