

3'-O-(5-Chloroanthraniloyl)-2'-deoxyadenosine-5'-O-triphosphate sodium

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: **Keep away from direct sunlight**
Store at -20°C
Actual storage temperature shall be subject to the COA.

Biological Description

Description	3'-O-(5-Chloroanthraniloyl)-2'-deoxyadenosine-5'-O-triphosphate (sodium) is a fluorescently labeled derivative of 2'-deoxyadenosine triphosphate (dATP) and can be used as a probe.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481