


1-(2-Aminoethyl)piperidin-4-ol

Chemical Properties

CAS No. :	129999-60-6
Formula:	C7H16N2O
Molecular Weight:	144.22
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	1-(2-Aminoethyl)piperidin-4-ol is a PROTAC linker utilized in the synthesis of PROTAC molecules.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	6.9339 mL	34.6693 mL	69.3385 mL
5 mM	1.3868 mL	6.9339 mL	13.8677 mL
10 mM	0.6934 mL	3.4669 mL	6.9339 mL
50 mM	0.1387 mL	0.6934 mL	1.3868 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481