

CEI12

Chemical Properties

CAS No. :	26046-92-4
Formula:	C72H97N13O15
Molecular Weight:	1384.64
Storage:	Keep away from moisture Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

Biological Description

Description	CEI12 is an angiotensin I-converting enzyme (ACEI) inhibitor peptide derived from α 1-casein peptides, suitable for hypertension research.
Targets(IC50)	Angiotensin-converting Enzyme (ACE)

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.7222 mL	3.611 mL	7.2221 mL
5 mM	0.1444 mL	0.7222 mL	1.4444 mL
10 mM	0.0722 mL	0.3611 mL	0.7222 mL
50 mM	0.0144 mL	0.0722 mL	0.1444 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481