

Macropa-PEG4-TPF

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

| | |
|---------------|--|
| Description | Macropa-PEG4-TPF (Compound 7g) is a bifunctional chelator. It binds to lysine residues on YS5, with an average of approximately 2.6 chelators per antibody. Macropa-PEG4-TPF is suitable for use in PET imaging studies. |
| Targets(IC50) | Others |

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481