

Rhodamine B-Dextran chloride (MW 200000)

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage:   
Keep away from direct sunlight  
Store at -20°C  
Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Rhodamine B-Dextran (MW 200000) is a dextran labeled with Rhodamine B, having a molecular weight of 200000. Rhodamine B-Chitosan (MW 200000) is significantly utilized in areas such as fluorescence labeling, cell imaging, vascular imaging, and drug delivery research with its excitation/emission at 546/610 nm.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481