

mmu-miR-139-3p Antagomir

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	mmu-miR-139-3p Antagomir is modified with a high-affinity cholesterol group. This miRNA antagomir inhibits the function of mature miRNA by specifically binding to it, thereby preventing its complementary pairing with target genes. Compared to miRNA inhibitors, miRNA antagomirs demonstrate greater stability and inhibitory effectiveness in animal studies, and they more easily penetrate cellular membranes and tissue gaps, concentrating in target cells.
Targets(IC50)	MicroRNA

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481