

10-Iododecan-1-ol

Chemical Properties

CAS No. :	57395-49-0
Formula:	C10H21IO
Molecular Weight:	284.18
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

Biological Description

Description	10-Iododecan-1-ol is a biochemical reagent.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.5189 mL	17.5945 mL	35.189 mL
5 mM	0.7038 mL	3.5189 mL	7.0378 mL
10 mM	0.3519 mL	1.7594 mL	3.5189 mL
50 mM	0.0704 mL	0.3519 mL	0.7038 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481