

## Me-Tetrazine-Me-Ph-amide-Lys(Boc)-Lys-N3

## Chemical Properties

CAS No. :	2861998-17-4
Formula:	C31H45N11O7
Molecular Weight:	683.77
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

## Biological Description

Description	Me-Tetrazine-Me-Ph-amide-Lys(Boc)-Lys-N3 is an ADC linker that facilitates the synthesis of antibody-drug conjugates (ADCs).
Targets(IC50)	ADC Linker

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.4625 mL	7.3124 mL	14.6248 mL
5 mM	0.2925 mL	1.4625 mL	2.925 mL
10 mM	0.1462 mL	0.7312 mL	1.4625 mL
50 mM	0.0292 mL	0.1462 mL	0.2925 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481