

## Cys(Npys)-(Arg)9 acetate

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Keep away from moisture

Storage:

Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Cys(Npys)-(Arg)9 acetate is a cell-penetrating peptide composed of 9 D-arginine residues and an activated cysteine residue C(Npys). This peptide acts as a carrier and is useful for studies involving conjugation and cellular permeability.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481