

TG (16:0_16:1_20:4)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

| | |
|---------------|-------------------------------------------------------------------------------------------------------|
| Description | TG (16:0_16:1_20:4) is a monoarachidic triglyceride, which is categorized as a triacylglycerol lipid. |
| Targets(IC50) | Endogenous Metabolite |

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481