


2,2'-(Piperazine-1,4-diyl)diethanol

Chemical Properties

CAS No. :	122-96-3
Formula:	C ₈ H ₁₈ N ₂ O ₂
Molecular Weight:	174.24
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	2,2'-(Piperazine-1,4-diyl)diethanol is a PROTAC linker used for synthesizing PROTAC molecules.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.7392 mL	28.6961 mL	57.3921 mL
5 mM	1.1478 mL	5.7392 mL	11.4784 mL
10 mM	0.5739 mL	2.8696 mL	5.7392 mL
50 mM	0.1148 mL	0.5739 mL	1.1478 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481