

YH-APN

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Keep away from direct sunlight
Store at -20°C
Actual storage temperature shall be subject to the COA.

Biological Description

Description	YH-APN is a high-performance APN-activatable near-infrared fluorescent probe. It offers high sensitivity, high selectivity, and rapid response. When administered through in situ spraying, it achieves an exceptional tumor-to-normal tissue signal-to-noise ratio. YH-APN specifically differentiates cancer cells with high APN expression from normal cells and can monitor cancer cell migration.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481