

8-iso-17-Phenyl trinor prostaglandin F2 α

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	8-Iso-17-phenyl trinor prostaglandin F2 α (8-Iso-17-phenyl trinor PGF2 α) is the C-8 epimer of Bimatoprost, functioning as an analog of PGF2 α . Bimatoprost is utilized in studies related to glaucoma and ocular hypertension.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481