

## Tisocalcitate

## Chemical Properties

CAS No. :	156965-06-9
Formula:	C <sub>31</sub> H <sub>48</sub> O <sub>5</sub>
Molecular Weight:	500.72
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

## Biological Description

Description	Tisocalcitate is a Vitamin D (Vitamin D) analog utilized in the study of plaque psoriasis.
-------------	--

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.9971 mL	9.9856 mL	19.9712 mL
5 mM	0.3994 mL	1.9971 mL	3.9942 mL
10 mM	0.1997 mL	0.9986 mL	1.9971 mL
50 mM	0.0399 mL	0.1997 mL	0.3994 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481