

## N3-PEG8-Oxydiacetamide-Phe-Lys(MMT)-PAB-PNP

## Chemical Properties

CAS No. :	1039045-00-5
Formula:	C71H89N9O19
Molecular Weight:	1372.54
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

## Biological Description

Description	N3-PEG8-Oxydiacetamide-Phe-Lys(MMT)-PAB-PNP [compound 18] is an intermediate used in the synthesis of SN-38 derivatives.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.7286 mL	3.6429 mL	7.2858 mL
5 mM	0.1457 mL	0.7286 mL	1.4572 mL
10 mM	0.0729 mL	0.3643 mL	0.7286 mL
50 mM	0.0146 mL	0.0729 mL	0.1457 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481