

Gadofosveset trisodium

Chemical Properties

CAS No. :	211570-55-7
Formula:	C33H40GdN3Na3O15P
Molecular Weight:	975.87
Storage:	Keep away from direct sunlight Store at -20°C <small>Actual storage temperature shall be subject to the COA.</small>

Biological Description

Description	Gadofosveset trisodium (MS 325) is an innovative contrast agent used in magnetic resonance angiography. It facilitates comprehensive high-resolution imaging of both arteries and veins. Gadofosveset trisodium is utilized for assessing venous, dynamic, and functional vascular disorders.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.0247 mL	5.1236 mL	10.2473 mL
5 mM	0.2049 mL	1.0247 mL	2.0495 mL
10 mM	0.1025 mL	0.5124 mL	1.0247 mL
50 mM	0.0205 mL	0.1025 mL	0.2049 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481