

## CY3-Alginate

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage:

Keep away from direct sunlight

Store at -20°C

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	CY3-Alginate is a sodium alginate labeled with CY3. It merges the polysaccharide properties of sodium alginate with the fluorescent features of CY3, making it significant in fluorescence labeling, cell imaging, drug delivery studies, and other areas (Ex/Em = 550/570 nm).
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481