

Rp-dCTPαS

Chemical Properties

CAS No. : 80951-76-4
Formula: C₉H₁₆N₃O₁₂P₃S
Molecular Weight: 483.22
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	Rp-dCTPαS is a monomeric material used in nucleic acid synthesis.
Targets(IC50)	Nucleoside Antimetabolite/Analog

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.0695 mL	10.3473 mL	20.6945 mL
5 mM	0.4139 mL	2.0695 mL	4.1389 mL
10 mM	0.2069 mL	1.0347 mL	2.0695 mL
50 mM	0.0414 mL	0.2069 mL	0.4139 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481