

Sp-2'-O-Me-GTP α S tetrasodium

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Sp-2'-O-Me-GTP α S tetrasodium is a GTP analog. It is utilized in enzyme-catalyzed processes to incorporate phosphorothioate and modified 2'-OH groups into RNA.
Targets(IC50)	Nucleoside Antimetabolite/Analog,DNA/RNA Synthesis

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481