

Alloc-NH-PEG8-PA

Chemical Properties

CAS No. :	2897625-33-9
Formula:	C ₂₃ H ₄₃ N ₁₂ O ₁₂
Molecular Weight:	525.59
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

Biological Description

Description	Alloc-NH-PEG8-PA (Compound 48) is a linker used for the conjugation of Pimodivir.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.9026 mL	9.5131 mL	19.0262 mL
5 mM	0.3805 mL	1.9026 mL	3.8052 mL
10 mM	0.1903 mL	0.9513 mL	1.9026 mL
50 mM	0.0381 mL	0.1903 mL	0.3805 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481