


Carbobenzoxy-4-aminobutyric acid

Chemical Properties

CAS No. :	5105-78-2
Formula:	C ₁₂ H ₁₅ NO ₄
Molecular Weight:	237.26
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	Carbobenzoxy-4-aminobutyric acid is a PROTAC linker used in the synthesis of PROTAC molecules.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.2148 mL	21.0739 mL	42.1479 mL
5 mM	0.843 mL	4.2148 mL	8.4296 mL
10 mM	0.4215 mL	2.1074 mL	4.2148 mL
50 mM	0.0843 mL	0.4215 mL	0.843 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481