

Human Serum Albumin-CY5.5

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: **Keep away from direct sunlight**
Store at -20°C
Actual storage temperature shall be subject to the COA.

Biological Description

Description	Human Serum Albumin-CY5.5 (HSA-CY5.5) is a CY5.5-labeled serum albumin. It plays a significant role in fields such as fluorescence labeling and drug delivery research, with excitation/emission wavelengths of 680/710 nm.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481