

## N-Succinyl-D-alanine

## Chemical Properties

CAS No. :	112243-89-7
Formula:	C7H11NO5
Molecular Weight:	189.17
Storage:	Keep away from moisture Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

## Biological Description

Description	N-Succinyl-D-alanine is a peptide ligand of the antibiotic Ristocetin A.
Targets(IC50)	Antibiotic

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.2863 mL	26.4313 mL	52.8625 mL
5 mM	1.0573 mL	5.2863 mL	10.5725 mL
10 mM	0.5286 mL	2.6431 mL	5.2863 mL
50 mM	0.1057 mL	0.5286 mL	1.0573 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481