

C6 Biotin ganglioside gm3 (d18:1/6:0) ammonium

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: **Keep away from direct sunlight**
Store at -20°C
Actual storage temperature shall be subject to the COA.

Biological Description

Description	C6 Biotin ganglioside gm3 (d18:1/6:0) (Biotin-C6:0-gm3) ammonium is a derivative of ganglioside GM3. The compound can be used for the detection or immobilization of ganglioside GM3 through interaction with biotin ligands.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481