

mmu-miR-139-3p Agomir

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	mmu-miR-139-3p Agomir is a specially chemically modified miRNA mimic, featuring methoxy modifications across the full length of the mature miRNA chain, with 2 and 4 thiophosphate linkages at the 5' and 3' ends, respectively, and a high-affinity cholesterol modification attached at the 3' end. This miRNA agomir mimics endogenous miRNA and enhances miRNA activity. Compared to regular mimics, miRNA agomirs show increased stability and inhibitory effects in animal experiments, and they more effectively penetrate cell membranes and tissue gaps to accumulate in target cells.
Targets(IC50)	MicroRNA

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481