

MC-Phe-Lys(MMT)-PAB-PNP

Chemical Properties

CAS No. :	159857-66-6
Formula:	C59H60N6O11
Molecular Weight:	1029.16
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

Biological Description

Description	MC-Phe-Lys(MMT)-PAB-PNP (Compound 27) is an intermediate that can be utilized in research related to cancer, autoimmune diseases, or infectious diseases.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.9717 mL	4.8583 mL	9.7167 mL
5 mM	0.1943 mL	0.9717 mL	1.9433 mL
10 mM	0.0972 mL	0.4858 mL	0.9717 mL
50 mM	0.0194 mL	0.0972 mL	0.1943 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481