

1,2-Dideoxy-D-ribofuranose

Chemical Properties

CAS No. :	91547-59-0
Formula:	C ₅ H ₁₀ O ₃
Molecular Weight:	118.131
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

Biological Description

Description	1,2-Dideoxy-D-ribofuranose is a monomeric material utilized in the synthesis of nucleic acids.
Targets(IC50)	Nucleoside Antimetabolite/Analog

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	8.4653 mL	42.3263 mL	84.6525 mL
5 mM	1.6931 mL	8.4653 mL	16.9305 mL
10 mM	0.8465 mL	4.2326 mL	8.4653 mL
50 mM	0.1693 mL	0.8465 mL	1.6931 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481